

200 mM NaCl
0.25% Triton X-100

- Solubilize the recombinant protein in solubilization buffer: expect \pm 10 mg/L culture and we want \pm 2 mg/ml final concentration. Resuspend the pellet by vortexing and/or sonicating: rotate for at least 1 hr at RT.

Inclusion body solubilization buffer: 50 mM Tris/HCl pH 8.0
200 mM NaCl
8 M Urea
5 mM Imidazole

- Prepare a SDS-PAGE gel of the appropriate % to estimate protein yield (see below).

- Remove any non-solubilized material by centrifugation: 12,000 x g for 15 min at 4°C. Add $\frac{1}{3}$ volume of 50 mM Tris/HCl pH 8.0, 200 mM NaCl, 5 mM Imidazole to reduce the urea concentration of the protein extraction SN to 6M. Run some of the protein extract out on SDS-PAGE together with known amounts of BSA to estimate the protein concentration. The sample can be stored at -80°C until further processing.

III Protein Purification

Ni-NTA from Invitrogen comes as a 50% solution in 30% EtOH, and can bind 5-10 mg of protein per ml of resin: calculate how much Ni-slurry you need. The following protocol is for a 5 ml bed volume column (10 ml of Ni-slurry) / batch-wise purification: do not allow the resin to dry at any point!

- Equilibrate the nickel-sepharose: 3 x 25 ml Milli-Q washes
(each wash is for 5 min) 3 x 25 ml NiCl₂ (5 mg/ml) washes
3 x 25 ml denaturing binding buffer

Denaturing binding buffer: 50 mM Tris/HCl pH 8.0
200 mM NaCl
6 M Urea
5 mM Imidazole

- Incubate the nickel-sepharose with the 6 M urea protein extract for 1 hr at RT by rotary shaking (take a 100 μ l sample before start: **START**). Afterwards, let the nickel-sepharose settle by gravity and remove SN + store at - 80°C (take a 100 μ l sample: **FLOW THROUGH**).

- Wash column 3 x with 25 ml denaturing wash buffer; each wash is for 5 min (take a 100 μ l sample of each wash: **WASH I, II and III**).

Denaturing wash buffer: 50 mM Tris/HCl pH 8.0
200 mM NaCl
6 M Urea
20 mM Imidazole

- Elute bound protein with 25 ml of denaturing elution buffer (take a 100 μ l sample: **ELUATE**).

Denaturing elution buffer: 50 mM Tris/HCl pH 8.0
100 mM NaCl
6 M Urea
250 mM Imidazole

The same resin can be used three to five times before EDTA stripping and nickel recharging are necessary (see 10.11.17 of Current Protocols in Molecular Biology). If resin is to be stored for longer than 1 day, wash with 10 bed volumes 20% EtOH and add 2 bed volumes 20% EtOH (make sure column is sealed to prevent evaporation). Only one protein should be purified on any given column.

- Concentrate eluate with a 3 kD cut-off Amicon concentrator to a final volume of \pm 5 ml.

Dialyse the protein O/N at 4°C in 2L of Denaturing dialysis buffer (50 mM Tris/HCl pH 8.0, 100 mM NaCl, 6 M Urea) to get rid of the imidazole.

- Prepare SDS-PAGE gel and assess protein recovery by staining with Coomassie Blue (include BSA standard to determine final yield). Alternatively, you can also estimate from O.D.₂₈₀ = 1, representing 1 mg/ml solution.

IV Protein refolding

I will assume a final yield of 30 mg protein (3L of culture) in 5 ml = 6 mg/ml!

- Rapidly dilute by drop wise addition to a final protein concentration of 100 μ g/ml in 295 ml Refolding buffer I:

50 mM Tris/HCl pH 8.0
100 mM NaCl
1 M Urea
5 mM CaCl₂
1 μ M ZnCl₂

Add the purified protein slowly from a pipet while stirring the refolding buffer to ensure thorough mixing and continue to stir at RT for 2 hrs. Store the refolding mixture at 4°C O/N without stirring.

- Concentrate refolding mixture with Vivacell 70 ml concentrators (5 kD cut-off) to \pm 5ml: afterwards clean concentrators + store in 20% EtOH so that they can be re-used. Rapidly dilute concentrate by drop wise addition to a final protein concentration of 100 μ g/ml in 295 ml Refolding buffer II:

50 mM Tris/HCl pH 8.0
100 mM NaCl
5 mM CaCl₂
1 μM ZnCl₂

Add the purified protein slowly from a pipet while stirring the refolding buffer to ensure thorough mixing and continue to stir at RT for 2 hrs. Store the refolding mixture at 4°C O/N without stirring.

- Concentrate refolding mixture with Vivacell 70 ml concentrators (5 kD cut-off) to ± 5ml. Dialyse the refolding mixture O/N at 4°C against 2L of Refolding buffer II to completely get rid of the urea.

- Perform a protein assay to determine the final protein concentration. Perform a α-casein zymogram and an enzyme assay to determine the specific activity (also run some out on a regular SDS-PAGE). Freeze 100 μl aliquots at -80°C.

! THE ABOVE REFOLDING PROTOCOL LEADS TO AUTO-DEGRADATION OF THE FULL LENGTH (CLEAVES OF ITS PRO- AND HEMOPEXIN DOMAINS) AND THE CATALYTIC DOMAIN (CLEAVES OF ITS PRO-DOMAIN): REPEAT IN THE PRESENCE OF EITHER A GENERAL MMP INHIBITOR SUCH AS 10 μM 1,10-PHENANTHROLINE OR AN INHIBITING ZINC CONCENTRATION OF 5 mM ZnCl₂. REMOVE INHIBITORS BY DIALYSIS ONCE REFOLDING IS COMPLETE!